

Ansell Touch n` Tuff Green Nitrile Gloves - 92-500 92-600 92-605

PROVEN SPLASH RESISTANCE AGAINST HAZARDOUS CHEMICALS

- Further testing of the TouchNTuff® by a certified body on an even wider range of chemicals confirms that it resists a greater variety of industrial chemicals for longer periods than any other nitrile disposable gloves.
- Manufactured from nitrile, TouchNTuff® offers up to four times the puncture resistance of comparable natural-latex gloves, and three times the resistance of similar neoprene gloves.
- It contains no natural rubber proteins, which means no risk of Type I allergies for the wearer. Primary skin irritation studies and Insult Patch tests have also shown no evidence of risk from irritation or allergic contact dermatitis.
- With a unique “Thin Nitrile Technology” formulation, this glove offers easy donning and strong grip in wet or dry conditions.
- The glove is highly versatile and suitable for many different uses.

REMARKS

- 92-600: Fit for special purpose
- TouchNTuff® is anti-staticity tested (EN 1149-1 & 1149-3) and meets the requirements defined in EN 1149-5

PRIMARY INDUSTRIES

- Agriculture & Viticulture
- Automotive & Transportation
- Building & Construction
- Chemical
- Electronics
- Food Processing & Service
- Glass
- Maintenance
- Metal Fabrication
- Pharmaceuticals

IDEAL APPLICATIONS


- Chemical handling
- Laboratory analysis
- Pharmaceuticals
- Paint and spray shops
- Printing industry
- Electronics
- Light assembly tasks
- Light assembly of oil-coated pieces
- Glass manufacturing
- Handling of cytostatics
- Intricate parts handling

SPECIFICATIONS

- Coating material: Nitrile
- Liner material: n/a
- Grip Design: Smooth
- Cuff Style: Rolled beaded

- Colour: Green
- AQL EN374: 1.5
- Powder-free: Yes
- EN Size: 6½-7, 7½-8, 8½-9, 9½-10
- Length mm: 240
- Thickness mm: 0.12

CATEGORY III



PERMEATION BREAKTHROUGH TIMES ACCORDING TO EN374-3:2003

Chemical Agent	CAS Number	Breakthrough Time	Protection Index	Notified Body	EN Standard
1-Methoxy-2-Propanol	107-98-2	14	1	Centexbel	374-3:2003
Isopropanol 70% (Ipasept)	67-63-0	178	4	Centexbel	374-3:2003
Acetic Acid, Glacial	64-19-7	7	0	Satra	374-3:2003
Acetonitrile 73% + Methyl Alcohol 25% + Ammonia 2%		1	0	Centexbel	374-3:2003
Acrylamide, 40%	79-06-1	> 480	6	Force Technology	374-3:2003
Allylchloride	107-05-1	< 5	0	Centexbel	374-3:2003
Ammonium Hydroxide, 25%	1336-21-6	29	1	Centexbel	374-3:2003
Anioxyde™ 1000	79-21-0	> 480	6	Force Technology	374-3:2003
Benzyl Alcohol	100-51-6	10	1	Centexbel	374-3:2003
Bromochloromethane	74-97-5	88	3	Centexbel	374-3:2003
Butyl Alcohol	71-36-3	56	2	Centexbel	374-3:2003
Cacodylic acid Sodium salt buffer 0,1M		> 480	6	Centexbel	374-3:2003
Caffeine 1.6%	58-08-2	> 480	6	Centexbel	374-3:2003
Carbon disulfide	75-15-0	< 5	0	Centexbel	374-3:2003
Chlorobutane	25154-42-1	< 5	0	Centexbel	374-3:2003
Chloroform	67-66-3	0	0	Centexbel	374-3:2003
Cidex™	111-30-8	> 480	6	Force Technology	374-3:2003
Cidex™ OPA	643-79-8	> 480	6	Force Technology	374-3:2003
Cyclohexane	110-82-7	> 480	6	Centexbel	374-3:2003
Dibromomethane	74-95-3	< 5	0	Centexbel	374-3:2003
Diesel fuel	68334-30-5	> 480	6	Centexbel	374-3:2003
Diethylamine	109-89-7	1	0	Centexbel	374-3:2003
Dimethyl Sulfoxide	67-68-5	5	0	Centexbel	374-3:2003
Dimethylformamide	68-12-2	< 5	0	Centexbel	374-3:2003
Ethanol, 70%	64-17-5	27	1	Centexbel	374-3:2003
Ethidium bromide in water (saturated, ± 5%)	1239-45-8	> 480	6	Centexbel	374-3:2003
Ethyl acetate 86% + Methyl Alcohol 9% + Ammonia 5%		1	0	Centexbel	374-3:2003
Formaldehyde 4% in Phosphatebuffer	50-00-0	> 480	6	Centexbel	374-3:2003
Formaldehyde, 35%	50-00-0	> 480	6	Centexbel	374-3:2003
Gasoline	8006-61-9	84	3	Centexbel	374-3:2003
Glutaraldehyde, 50%	111-30-8	> 480	6	Centexbel	374-3:2003
Glutaric dialdehyde 2,5%, cacodylic acid, sodium salt		> 480	6	Centexbel	374-3:2003
Heptane 98% + 1-butyl alcohol 2%	142-82-5	9	0	Centexbel	374-3:2003
Heptane 98% + 3-methyl-1-butyl alcohol 2%	142-82-5	16	1	Centexbel	374-3:2003
Hexane	110-54-3	> 480	6	Centexbel	374-3:2003
Hydrofluoric Acid, 48%	7664-39-3	< 5	0	Centexbel	374-3:2003

█ Hydrogen Bromide, 49%	10035-10-6	> 480	6	Centexbel	374-3:2003
█ Hydrogen Peroxide, 30%	7722-84-1	41	2	Centexbel	374-3:2003
█ Iso-Octane	540-84-1	> 480	6	Centexbel	374-3:2003
█ Kerosene	64742-81-0	> 480	6	Centexbel	374-3:2003
█ Methanol	67-56-1	1	0	Centexbel	374-3:2003
█ Methyl ethyl ketone	78-93-3	< 5	0	Centexbel	374-3:2003
█ Methyl Isobutyl Ketone	108-10-1	1	0	Centexbel	374-3:2003
█ Methyl Sulfoxide 5% in Citratebuffer		> 480	6	Centexbel	374-3:2003
█ Methylviolet 1%	8004-87-3	> 480	6	Centexbel	374-3:2003
█ n-Undecane	1120-21-4	> 480	6	Centexbel	374-3:2003
█ Nitric Acid, 70%	7697-37-2	< 5	0	Centexbel	374-3:2003
█ Peracetic acid, 39%	79-21-0	9	0	Force Technology	374-3:2003
█ Perchloroethylene	127-18-4	8	0	Centexbel	374-3:2003
█ Potassium permanganate 5%	7722-64-7	120	4	Centexbel	374-3:2003
█ Salicylic acid 2% in Peanut oil		> 480	6	Centexbel	374-3:2003
█ Sodium Hydroxide, 50%	1310-73-2	> 480	6	Centexbel	374-3:2003
█ Sulphuric acid, 50%	7664-93-9	> 480	6	Centexbel	374-3:2003
█ Sulphuric acid, 99-100%	7664-93-9	1	0	Centexbel	374-3:2003
█ Tetrahydrofuran	109-99-9	< 5	0	Centexbel	374-3:2003
█ Toluene	108-88-3	1	0	Centexbel	374-3:2003
█ Triethylamine	121-44-8	155	4	Centexbel	374-3:2003
█ 1,1,1-trichloro-2-methyl-2-propyl alcohol in Peanut oil		> 480	6	Centexbel	374-3:2003
█ Formaldehyde, 24.5%	50-00-0	> 480	6	Centexbel	374-3:2003
█ Acrylic Acid	79-10-7	< 5	0	Centexbel	374-3:2003
█ Dichloroethane		< 1	0	Centexbel	374-3:2003
█ Cyclohexanone	108-94-1	< 5	0	Centexbel	374-3:2003
█ Dibromoethane	106-93-4	< 1	0	Centexbel	374-3:2003
█ Diethyl ether	60-29-7	< 1	0	Centexbel	374-3:2003
█ Ditranol 0,7% in liquid paraffin thin		1.6	0	Centexbel	374-3:2003
█ Ethanol, 95%	64-17-5	16	1	Centexbel	374-3:2003
█ Ethyl Acetate	141-78-6	1	0	Centexbel	374-3:2003
█ Heptane	142-82-5	> 480	6	Centexbel	374-3:2003
█ Hydrochloric Acid, 37%	7647-01-0	51	2	Centexbel	374-3:2003
█ Isopropanol	67-63-0	117	3	Centexbel	374-3:2003
█ Methyl sulfoxide 20% in RPMI 1640 culture 80%		> 480	6	Centexbel	374-3:2003
█ Methyl-t-butyl Ether	1634-04-4	14	1	Centexbel	374-3:2003
█ Methylmethacrylate	80-62-6	2	0	Force Technology	374-3:2003
█ Nicotine	54-11-5	25	1	Force Technology	374-3:2003
█ Nitric Acid, 50%	7697-37-2	9	0	Centexbel	374-3:2003
█ Tetrahydrofuran/n-Heptan, ratio:60%-40%		<5	0	Centexbel	374-3:2003
█ White Spirit	64742-88-7	285	5	Centexbel	374-3:2003
█ Xylene	1330-20-7	< 5	0	Centexbel	374-3:2003