

Terumo® 3-part Syringe

Product Information Card



Description:

Sterile 3-part syringe for single use of *Terumo (Philippines) Corporation origin*.

The syringe is intended for the injection of fluids - immediately after filling - into the human body, for aspiration of fluids or for the sampling of blood or other body fluids.

Components:

Description	Raw Materials
Barrel	Polypropylene
Plunger	Polystyrene (for 1 ml syringe only) Polypropylene
Gasket	Thermoplastic elastomer
Lubricant	Silicone oil



Product Range:

Product Code	Size	Tip	Replacement for
SS+01T1	1 ml	Concentric, Luer Slip	BS-01T
SS+02S1	2.5 ml	Concentric, Luer Slip	SS*02SE1
SS+03L1	3 ml	Concentric, Luer Lock	
SS+05S1	5 ml	Concentric, Luer Slip	
SS+10ES1	10 ml	Eccentric, Luer Slip	SS-10ES
SS+20ES1	20 ml	Eccentric, Luer Slip	
SS+20L1	20 ml	Concentric, Luer Lock	
SS+50ES1	50 ml	Eccentric, Luer Slip	BS-50ES
SS+50L1	50 ml	Concentric, Luer Lock	BS-50LG
SS+50C1	50 ml	Catheter Tip	BS-50C

Technical Data:

Product Code	Nominal Capacity	Graduated Capacity	Increments
SS+01T1	1 ml	1 ml	0.01
SS+02S1	2.5 ml	2.5 ml	0.1
SS+03L1	3 ml	3 ml	0.1
SS+05S1	5 ml	6 ml	0.2
SS+10ES1	10 ml	10 ml	0.2
SS+20ES1	20 ml	20 ml	1.0
SS+20L1	20 ml	20 ml	1.0
SS+50ES1	50 ml	60 ml	1.0
SS+50L1	50 ml	60 ml	1.0
SS+50C1	50 ml	60 ml	1.0

Terumo® 3-Part Syringe

Packaging:

Single unit: blister pack (film/film)

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SS+01T1	100	1800	7.50 kg	430 x 418 x 306 mm
SS+02S1	100	2400	7.00 kg	496 x 427 x 427 mm
SS+03L1	100	1800	7.50 kg	640 x 385 x 262 mm
SS+05S1	100	1200	9.40 kg	480 x 379 x 362 mm
SS+10ES1	100	1200	11.00 kg	546 x 454 x 428 mm
SS+20ES1	50	600	10.15 kg	570 x 400 x 366 mm
SS+20L1	50	600	10.15 kg	570 x 400 x 366 mm
SS+50ES1	25	100	5.10 kg	466 x 312 x 268 mm
SS+50L1	25	100	5.10 kg	466 x 312 x 268 mm
SS+50C1	25	100	5.40 kg	466 x 312 x 268 mm

Sterilization:

With e-beam radiation to a Sterility Assurance Level of at least 10⁻⁶

Storage:

Shelf life of 5 years

* The batch number and expiry date are embossed in the film of the blister pack

* Do not store at extreme temperature and humidity

Precautions for Use:

Check integrity of the unit package before use

For single use only

Dispose of safely after single use to avoid risk of infection

Status of Regulation:

Class IIa, sterile product following MDD 93/42/EEC

Notified body:

TÜV Rheinland 0197

Manufacturer:

TERUMO (PHILIPPINES) CORPORATION

EC representative:

TERUMO EUROPE NV

Interleuvenlaan 40

B - 3001 Leuven

Belgium



Terumo Corporation

44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072, Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

Terumo Europe NV

Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

Terumo Europe NV Benelux Sales Division

Belgium:
Tel.: +32 16 39 25 80
The Netherlands:
Tel.: 0800 0220396

Terumo Europe NV

Emerging Markets Division
Tel.: +32 16 38 13 08

Terumo Deutschland GmbH

Tel.: +49 6196 80 230

Terumo Deutschland GmbH

Zweigniederlassung
Spreitenbach
Tel.: +41 56 419 10 10

Laboratoires

Terumo France S.A.S.
Tel.: +33 130 96 13 00

Terumo Italia S.r.l.

Tel.: +39 0651 96 14 20

Terumo Europe España SL

Tel.: +34 902 10 12 98

Terumo UK Ltd

Tel.: +44 1276 480 440

Terumo Russia LLC

Tel.: +7 495 988 4740

Terumo Sweden AB

Tel.: +46 3174 85 880

Terumo Denmark

Tel.: +45 7020 93 80

Terumo Corporation Dubai Branch (Middle East and Africa)

Tel.: +971 4 221 2220

Terumo BCT Tıbbi Cihazlar Dağıtım ve Hizmetleri A.Ş.

Tel.: +90 216 590 0074



HO 20GB-0307TE-VIII(01.15)E

© Registered Trademark
Published by Terumo Europe NV