

Sani-Cloth® AF is an alcohol-free, quaternary ammonium formulation providing a two-in-one solution ideal for cleaning and disinfecting.

Tuberculocidal, virucidal, bactericidal, and yeasticidal in 2 minutes.



Recommended Uses:

- + Routine cleaning & disinfection
- + High-footfall areas
- + High touch point surfaces
- + Common equipment & non-invasive medical devices

+ Ready to Use

Pre-dosed, alcohol-free, low linting disposable cleaning and disinfection wipes, requiring no mixing or measuring.

+ Free From

Alcohol free, bleach free and PHMB free, causing no skin irritation or sensitisation.

+ Improved Efficiencies

One product to clean and disinfect, providing you with consistent and reliable results.

+ Supporting Resources

Point of Use accessories, including canister bracket and soft pack holders, "I Am Clean" tape & sticky notes, IFU videos and posters.

+ Proven Efficacy

Complete disinfection of all listed pathogens in 2 minutes (including norovirus in 30 seconds).

+ Cost Effective

Using the same solution for cleaning and disinfecting saves you time and money.

+ Maximum Flexibility

Suitable for any healthcare environment available in soft pack, canisters and buckets.

+ Supports Carbon Net Zero Goals

Manufactured in the UK, resulting in supply chain reliability and carbon footprint reduction.

Product Efficacy

| | Test Method | Organism | Contact Time | |
|---------------------|-------------|-----------------------------------|--------------|------------|
| Bacteria | EN 16615 | Pseudomonas aeruginosa | 60 sec (c) | 2 min (d) |
| | EN 16615 | Staphylococcus aureus | 60 sec (c) | 2 min (d) |
| | EN 16615 | Enterococcus hirae | 60 sec (c) | 2 min (d) |
| | EN 13727 | Pseudomonas aeruginosa | 30 sec (c) | |
| | EN 13727 | Staphylococcus aureus | 30 sec (c) | |
| | EN 13727 | Enterococcus hirae | 30 sec (c) | |
| Yeast | EN 16615 | Candida albicans | 60 sec (c) | 2 min (d) |
| | EN 13624 | Candida albicans | 30 sec (c) | |
| | EN 13624 | Candida auris | 30 sec (c) | |
| Enveloped virus | EN 14476 | Vaccinia Virus | 30 sec (c) | 60 sec (d) |
| | EN 14476 | Corona, HIV, Hep B & C, SARS, Flu | 30 sec (c) | 2 min (d) |
| | EN 14476 | Human Coronavirus ATCC VR-740 | 30 sec (c) | |
| Non-enveloped virus | EN 14476 | Norovirus | 30 sec (c) | 60 sec (d) |
| TB | EN 14348 | Mycobacterium terrae | 2 min (c) | |

*Clean Conditions (c) Dirty Conditions (d)



Product Details

| Format | Wipe Size | Wipe Count | QTY per Case | Classification | Product Code |
|-----------------|---------------|------------|--------------|----------------|--------------|
| 100ct Soft Pack | 280mm x 200mm | 100 | 12 | MD Class IIa | EAF1100 |
| 225ct Bucket | 270mm x 267mm | 225 | 1 | MD Class IIa | EAF1110 |
| 200ct Soft Pack | 267mm x 210mm | 200 | 6 | MD Class IIa | EAF1130 |
| 40ct Soft Pack | 200mm x 200mm | 40 | 24 | MD Class IIa | EAF1140 |



For more information, write to us: contact@pdi-emea.com

T: + 44 (0) 8081 697 945 E: contact@pdi-emea.com W: pdi-intl.com

PDI International