


MD300C19



FINGERTIP PULSE OXIMETER

Pulse oximeter is a very important and common device to check oxygen saturation (SpO₂) and pulse rate. As a small, compact, simple, reliable and durable physiological monitoring device, Fingertip pulse oximeter is widely applied in clinic, hospital, social medical organization, and homecare etc.



FEATURES

- ♥ High brightness LED displays SpO₂, PR, and Pulse bar
- ♥ Low power consumption; battery-low indication
- ♥ Automatically power off; 2pcs AAA-size batteries
- ♥ Two display modes
- ♥ Certifications: CE, FDA, CFDA



TECHNICAL SPECIFICATIONS

Application Range

Adult, adolescent and child

| SpO ₂ | | Pulse Rate | |
|--------------------------|---------------------------------------|-----------------------|---|
| Measurement Range: | 70%~100% | Measurement Range: | 30bpm~250bpm |
| Resolution: | 1% | Resolution: | 1bpm |
| Measurement Accuracy: | 70%~100% ±2%; 0%~69% no definition | Measurement Accuracy: | 30bpm~99bpm ±2bpm; 100bpm~250bpm ±2% |
| Display | | Dimensions and Weight | |
| Type: | LED | Dimension: | 58mm (L) * 32mm (W) * 34mm (H) |
| Parameters: | SpO ₂ , PR, Pulse bar | Weight: | 50g (including batteries) |
| Power Requirements | | | |
| Battery: | Two AAA-size 1.5V Batteries | Power Consumption: | Less than 40mA |
| Environment Requirements | | | |
| Operation Temperature: | 5°C~40°C | Operation Humidity: | 15% ~ 93%, no condensation |
| Storage Temperature: | -25°C~+70°C | Storage Humidity: | ≤93%, no condensation |

Beijing Choice Electronic Tech Co., Ltd.

Room 4104, No. A12 Yuquan Road, Haidian District,
Beijing 100143, P.R.China
Website: www.choicemed.com Fax: 86 10-88204632
Email: market@choicemed.com Tel: 86 10-88204630/31

ChoiceMMed
Providing a Better Life

NOTE The specifications are subject to changes without prior notice.
Revised Date: 25 January 2019